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Giltig från: 2025-11-17

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
Aktiv befuktning, vuxen

<p>NIV Värmaren i "NIV-läge" Ingen etCO₂</p> <p>Nebulisatorn vid värmaren.</p>	 A photograph showing a white nebulizer connected to a blue corrugated tube. The tube is connected to a white and blue NIV (Non-Invasive Ventilation) mask. The nebulizer has a red label with the word 'Nebulizer'.
<p>Tubad – öppet sugsystem Snorkel + knä EtCO₂ via expiratorisk (vit slang) och grön adapter. Nebulisatorn vid värmaren.</p> <p><i>Deadspace ca 62 ml</i></p>	 A photograph showing a nebulizer connected to a blue corrugated tube. The tube is connected to a green and white adapter. A white tube (expiratorisk slang) is connected to the adapter. A green adapter is also visible. A small white label with a barcode is attached to the blue tube.
<p>Tubad + slutet sugsystem EtCO₂ som ovan. Nebulisatorn vid värmaren.</p> <p><i>Deadspace 62 ml</i></p>	 A photograph showing a nebulizer connected to a blue corrugated tube. The tube is connected to a green and white adapter. A white tube (expiratorisk slang) is connected to the adapter. A green adapter is also visible. A small white label with a barcode is attached to the blue tube.
<p>HFNC vuxen Y-koppling + vit slang kopplas bort. HFNC-grimma på blå slang Nebulisatorn vid värmaren. Tänk på ev dubbel dos inhalation.</p>	 A photograph showing a nebulizer connected to a blue corrugated tube. The tube is connected to a white Y-connector. A white corrugated tube (HFNC-grimma) is connected to the Y-connector. A green adapter is also visible.

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


Hamilton transportventilator

<p>Tubad vuxen + HME</p> <p><i>Deadspace 180 ml</i></p>	 A photograph showing a Hamilton transport ventilator setup for an adult patient with a Heat and Moisture Exchanger (HME) filter. The setup includes a clear plastic ventilator body, a blue corrugated breathing tube, and a yellow HME filter. A blue connector is visible on the right side.
<p>NIV vuxen</p>	 A photograph showing a Hamilton transport ventilator setup for an adult patient using Non-Invasive Ventilation (NIV). The setup includes a clear plastic ventilator body, a blue corrugated breathing tube, and a blue NIV mask. A blue connector is visible on the right side.
<p>Tubad barn 10-30 kg</p> <p><i>Deadspace ca 49 ml</i></p>	 A photograph showing a Hamilton transport ventilator setup for a child weighing 10-30 kg. The setup includes a clear plastic ventilator body, a blue corrugated breathing tube, and a yellow HME filter. A blue connector is visible on the right side.
<p>Tubad barn < 10 kg (Vt<60ml)</p> <p>Ev med knä, men det ökar deadspace.</p> <p>Obs! Vuxen CO2-adapter är tyvärr det enda som passar till neo-HME-filter.</p> <p><i>Deadspace ca 10 ml</i></p>	 A photograph showing a Hamilton transport ventilator setup for a child weighing less than 10 kg. The setup includes a clear plastic ventilator body, a blue corrugated breathing tube, and a blue connector. A small, clear plastic adapter is shown above the main setup.

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3. Maquet Critical Care AB, *Bruksanvisning SERVO-U Ventilator System v4.5, Rev.18*, publicerad 26 januari 2023