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Referenser

1. <https://www.sfog.se/media/336641/sfog-raad-transducerskydd-och-desinfektion-av-ultraljudsprober-vid-all-gynekologisk-och-obstetrisk-ultraljudsundersokning-2019-09-27.pdf>
2. Jacques S. Abramowicz, J S, et al. Guidelines for Cleaning Transvaginal Ultrasound Transducers Between Patients. WFUMB Safety Committee. *Ultrasound in Medicine & Biology* (2017), 43(5), 1076-9.
3. Nyhsen, C M, et al. Infection prevention and control in ultrasound - best practice recommendations from the European Society of Radiology Ultrasound Working Group. *Insights Imaging* (2017), 8, 523.
4. McDonnelli D, Russell, AD. Antiseptics and Disinfectants: Activity, Action, and Resistance. *Clin Microbiol Reviews* (1999), 12(1), 147–79.
5. NHS Scotland. Guidance for Decontamination of Semi-Critical Ultrasound Probes; Semi-invasive and Non-invasive Ultrasound Probes (2016).
6. Government of Western Australia, Dep of Health. Operational Directive 0404/12: Prevention of Cross infection in Diagnostic Ultrasound (2012).
7. Leigh J, et al. The cost of being clean: A cost analysis of nasopharyngoscope reprocessing techniques. *Laryngoscope* (2018) 128(1), 64-71.
8. Hitchcock B, et al. A randomised, single-blind comparison of high-level disinfectants for flexible nasendoscopes. *The Journal of Laryngology & Otology* (2016), 130(11), 983-989.

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Bilaga 2

Tristel säkerhet

Tristel Duo for Surfaces Foam Activated Foam Safety Data Sheet

1. Identification

Product Name:	Tristel Duo for Surfaces
Product Code:	DUO/SUR
Type of product:	Sporicidal cleaner for use on medical equipment and surfaces
Manufacturer:	Tristel Solutions Limited Lynx Business Park Fordham Road Snailwell Cambs, CB8 7NY, United Kingdom
Telephone Number:	+44 (0) 1638 721500 +44 (0) 7798 805692 (Out of hours contact)
Fax Number:	+44 (0) 1638 721911
Email:	healthandsafety@tristel.com
Chemical type:	Mild oxidising solution

2. Hazards

identification:	Chlorine dioxide generator WEL - 0.84 mg/m ³ short term, 15 minute reference period - 0.28 mg/m ³ long term, 8 hour TWA reference period
-----------------	---

3. Composition:

Ingredients	CAS No	EINECS No	Wt/Vol %	Symbol
Chlorine dioxide in aqueous solution	10049-04-4	233-162-8	0.1-0.12	ClO ₂
Water				H ₂ O

4. First aid measures:

Inhalation:	Non-toxic
Eye contact:	Rinse eyes with water
Skin contact:	Wash affected area with water
Ingestion:	Do not induce vomiting. Give milk or water to drink Seek medical advice where necessary

5. Fire fighting measures:

Non flammable

6. Accidental release measures:

Environmental precautions:	Environmental precautions required but product is biodegradable under OECD conditions operational 6/1995
Clean up method:	Dispose of to clinical waste
Clothing for disposal:	Wear appropriate gloves and apron

This document conforms
with the Regulations of
REACH, 2007

Tristel Duo for Surfaces Foam Activated Foam Safety Data Sheet

7. Handling and storage:

Handling guidelines:	Safe handling in accordance with label instructions Not to be mixed with other chemicals Keep from children
Storage guidelines:	Store within a temperature range of 10° C to 35° C Single use solution

8. Exposure controls/personal protection:

Personal protection:	Avoid eye contact and prolonged skin contact Gloves and apron recommended
Skin contact:	Low risk
Eye contact:	Low risk
Inhalation:	Low risk
Ingestion:	Low risk, substantial ingestion may cause discomfort to mouth and digestive tissues

9. Physical and chemical properties:

Physical state:	Liquid (Foam)
Appearance and odour:	White foam with mild odour
Evaporation rate:	As water
Boiling point:	As water
Freezing point:	As water
% Volatile (by weight):	Not known
Solubility in water (20°C):	Soluble
pH:	5.0 approximately
Specific gravity:	Not applicable

This document conforms
with the Regulations of
REACH, 2007

Tristel Duo for Surfaces Foam
Foam (Activator) & Foam (Base) Safety Data Sheet

1. Identification Product Name: Tristel Duo for Surfaces – Activator Solution Chemical Type: Sodium salt solution	1. Identification Product Name: Tristel Duo for Surfaces – Base Solution Chemical Type: Organic acid blend
2. Hazards Identification: No specific hazards Contact with acids liberates chlorine dioxide	2. Hazards Identification: No specific hazards
3. Composition/information on ingredients: Major ingredients: 0.5% Sodium Chlorite solution 2.0% Foaming Surfactant in de-mineralised water	3. Composition/information on ingredients: Major ingredients: 5% solution of Citric Acid 0.05% Boric Acid 0.05% Sorbic Acid 2.0% Foaming surfactant in de-mineralised water
4. First-aid measures: Inhalation: Non-toxic Eye contact: Rinse eyes with water Skin contact: Wash affected area with water Ingestion: Do not induce vomiting Give milk or water to drink Seek medical advice where necessary	4. First-aid measures: Inhalation: Non-toxic Eye contact: Rinse eyes with water Skin contact: Wash affected area with water Ingestion: Do not induce vomiting Give milk or water to drink Seek medical advice where necessary
5. Fire fighting measures: Non flammable	5. Fire fighting measures: Non flammable
6. Accidental release measures: Environmental Precautions: Environmental precautions required but product is biodegradable under OECD conditions operational 6/1995 Clean up method: Dispose of to clinical waste Clothing for disposal: Wear appropriate gloves	6. Accidental release measures: Environmental Precautions: Environmental precautions required but product is biodegradable under OECD conditions operational 6/1995 Clean up method: Dispose of to clinical waste Clothing for disposal: Wear appropriate gloves <div style="text-align: right; margin-top: 20px;"> This document conforms with the Regulations of REACH, 2007 </div>

Tristel Duo for Surfaces Foam
Foam (Activator) & Foam (Base) Safety Data Sheet

Product Name: Tristel Duo for Surfaces – Activator Solution	Product Name: Tristel Duo for Surfaces – Base Solution
7. Handling and storage:	7. Handling and storage:
<p>Handling guidelines: Safe handling in accordance with label instructions</p> <p>Not to be mixed with other chemicals</p> <p>Keep from children</p> <p>Storage guidelines: Store within a temperature range of 10° C to 35° C</p> <p>Shelf life – two years – see product for expiry date</p>	<p>Handling guidelines: Safe handling in accordance with label instructions</p> <p>Not to be mixed with other chemicals</p> <p>Keep from children</p> <p>Storage guidelines: Store within a temperature range of 10° C to 35° C</p> <p>Shelf life – two years – see product for expiry date</p>
8. Exposure controls/personal protection:	8. Exposure controls/personal protection:
<p>Personal protection: Avoid eye contact and prolonged skin contact</p> <p>Gloves and apron recommended</p> <p>Skin contact: Low risk</p> <p>Eye contact: Low risk</p> <p>Inhalation: Low risk</p> <p>Ingestion: Low risk, substantial ingestion may cause discomfort to mouth and digestive tissues</p>	<p>Personal protection: Avoid eye contact and prolonged skin contact</p> <p>Gloves and apron recommended</p> <p>Skin contact: Low risk</p> <p>Eye contact: Low risk</p> <p>Inhalation: Low risk</p> <p>Ingestion: Low risk, substantial ingestion may cause discomfort to mouth and digestive tissues</p>
9. Physical and chemical Properties:	9. Physical and chemical Properties:
<p>Physical state: Liquid</p> <p>Appearance and odour: Clear, no odour</p> <p>Evaporation rate: As water</p> <p>Boiling point: As water</p> <p>Freezing point: As water</p> <p>% Volatile (by weight): Not known</p> <p>Solubility in water (20° C): Soluble</p> <p>pH: 12.0 approximately</p> <p>Specific gravity: 1.005 @ 20° C</p>	<p>Physical state: Liquid</p> <p>Appearance and odour: Blue, no odour</p> <p>Evaporation rate: As water</p> <p>Boiling point: As water</p> <p>Freezing point: As water</p> <p>% Volatile (by weight): Not known</p> <p>Solubility in water (20° C): Soluble</p> <p>pH: 2.5 approximately</p> <p>Specific gravity: 1.025 @ 20° C</p>
	This document conforms with the Regulations of REACH, 2007

**Tristel Duo for Surfaces Foam
Foam (Activator) & Foam (Base) Safety Data Sheet**

Product Name: Tristel Duo for Surfaces – Activator Solution	Product Name: Tristel Duo for Surfaces – Base Solution
10. Stability and reactivity: No decomposition if stored and used as directed Hazardous decomposition products: None under normal use Chlorine donors and oxygen produced if heated	10. Stability and reactivity: No decomposition if stored and used as directed Hazardous decomposition products: None under normal use Not compatible with alkaline substance and chlorine donors
11. Toxicological information: Sodium Chlorite Acute Oral LD50: >4000 mg/kg Irritation to skin: Minimally irritating Irritation to eyes: Mild irritant Surfactant Acute Oral LD50: >2000 mg/kg Irritation to skin: Mildly irritating Irritation to eyes: Irritant	11. Toxicological information: Citric Acid Acute oral LD50: >11750 mg/kg Irritation to Skin: Mild irritant under prolonged contact Irritation to Eyes: Irritant Surfactant Acute Oral LD50: >2000 mg/kg Irritation to skin: Mildly irritating Irritation to eyes: Irritant
12. Ecological information: No known adverse effects from normal use	12. Ecological information: No known adverse effects from normal use
13. Disposal considerations: Packaging: Can be disposed of as normal waste in accordance with local authority regulations Contaminated packaging: May be disposed of safely under normal conditions in accordance with local authority regulations Product: Solution to be disposed of in accordance with spillage instructions as explained in accidental release measures	13. Disposal considerations: Packaging: Can be disposed of as normal waste in accordance with local authority regulations Contaminated packaging: May be disposed of safely under normal conditions in accordance with local authority regulations Product: Solution to be disposed of in accordance with spillage instructions as explained in accidental release measures
14. Transport information: Not hazardous for transport	14. Transport information: Not hazardous for transport <p style="text-align: right;">This document conforms with the Regulations of REACH, 2007</p>