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Referenser

1. <https://www.sfog.se/media/336641/sfog-raad-transducerskydd-och-desinfektion-av-ultraljudsprober-vid-all-gynekologisk-och-obstetrisk-ultraljudsundersokning-2019-09-27.pdf>
2. Jacques S. Abramowicz, J S, et al. Guidelines for Cleaning Transvaginal Ultrasound Transducers Between Patients. WFUMB Safety Committee. *Ultrasound in Medicine & Biology* (2017), 43(5), 1076-9.
3. Nyhsen, C M, et al. Infection prevention and control in ultrasound - best practice recommendations from the European Society of Radiology Ultrasound Working Group. *Insights Imaging* (2017), 8, 523.
4. McDonnelli D, Russell, AD. Antiseptics and Disinfectants: Activity, Action, and Resistance. *Clin Microbiol Reviews* (1999), 12(1), 147–79.
5. NHS Scotland. Guidance for Decontamination of Semi-Critical Ultrasound Probes; Semi-invasive and Non-invasive Ultrasound Probes (2016).
6. Government of Western Australia, Dep of Health. Operational Directive 0404/12: Prevention of Cross infection in Diagnostic Ultrasound (2012).
7. Leigh J, et al. The cost of being clean: A cost analysis of nasopharyngoscope reprocessing techniques. *Laryngoscope* (2018) 128(1), 64-71.
8. Hitchcock B, et al. A randomised, single-blind comparison of high-level disinfectants for flexible nasendoscopes. *The Journal of Laryngology & Otology* (2016), 130(11), 983-989.

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Bilaga 2

Tristel säkerhet

Tristel Duo for Surfaces Foam Activated Foam Safety Data Sheet

1. Identification

Product Name:	Tristel Duo for Surfaces
Product Code:	DUO/SUR
Type of product:	Sporicidal cleaner for use on medical equipment and surfaces
Manufacturer:	Tristel Solutions Limited Lynx Business Park Fordham Road Snailwell Cambs, CB8 7NY, United Kingdom
Telephone Number:	+44 (0) 1638 721500 +44 (0) 7798 805692 (Out of hours contact)
Fax Number:	+44 (0) 1638 721911
Email:	healthandsafety@tristel.com
Chemical type:	Mild oxidising solution

2. Hazards

identification:	Chlorine dioxide generator WEL - 0.84 mg/m ³ short term, 15 minute reference period - 0.28 mg/m ³ long term, 8 hour TWA reference period
-----------------	---

3. Composition:

Ingredients	CAS No	EINECS No	Wt/Vol %	Symbol
Chlorine dioxide in aqueous solution	10049-04-4	233-162-8	0.1-0.12	ClO ₂
Water				H ₂ O

4. First aid measures:

Inhalation:	Non-toxic
Eye contact:	Rinse eyes with water
Skin contact:	Wash affected area with water
Ingestion:	Do not induce vomiting. Give milk or water to drink Seek medical advice where necessary

5. Fire fighting measures:

Non flammable

6. Accidental release measures:

Environmental precautions:	Environmental precautions required but product is biodegradable under OECD conditions operational 6/1995
Clean up method:	Dispose of to clinical waste
Clothing for disposal:	Wear appropriate gloves and apron

This document conforms
with the Regulations of
REACH, 2007

Tristel Duo for Surfaces Foam Activated Foam Safety Data Sheet

7. Handling and storage:

Handling guidelines:	Safe handling in accordance with label instructions Not to be mixed with other chemicals Keep from children
Storage guidelines:	Store within a temperature range of 10° C to 35° C Single use solution

8. Exposure controls/personal protection:

Personal protection:	Avoid eye contact and prolonged skin contact Gloves and apron recommended
Skin contact:	Low risk
Eye contact:	Low risk
Inhalation:	Low risk
Ingestion:	Low risk, substantial ingestion may cause discomfort to mouth and digestive tissues

9. Physical and chemical properties:

Physical state:	Liquid (Foam)
Appearance and odour:	White foam with mild odour
Evaporation rate:	As water
Boiling point:	As water
Freezing point:	As water
% Volatile (by weight):	Not known
Solubility in water (20°C):	Soluble
pH:	5.0 approximately
Specific gravity:	Not applicable

This document conforms
with the Regulations of
REACH, 2007

Tristel Duo for Surfaces Foam Activated Foam Safety Data Sheet

10. Stability and reactivity:

Solution decomposes to simple salt solution

Hazardous decomposition products: None under normal use

Chlorine donors and oxygen produced if heated

11. Toxicological information:

Acute Oral LD ₅₀	>5000 mg/kg
Irritation to skin	Negative
Irritation to eyes	Negative
Sensitisation	Negative

12. Ecological information:

Presents no known hazards to the environment

13. Disposal considerations:

Packaging:	Can be disposed of as normal waste in accordance with local authority regulations
Contaminated packaging:	May be disposed of safely under normal conditions in accordance with local authority regulations
Product:	Dispose of to clinical waste

14. Transport information:

Not applicable

15. Regulatory information:

Not a licensed medicine

Approved under EEC 93/42 for contact with medical devices

16. Other information:

Safety phrases:

(2) Keep out of reach of children

(24 / 25) Avoid contact with skin and eyes

(50) Do not mix with other chemicals

This document conforms
with the Regulations of
REACH, 2007

Tristel Duo for Surfaces Foam Foam (Activator) & Foam (Base) Safety Data Sheet

1. Identification Product Name: Tristel Duo for Surfaces – Activator Solution Chemical Type: Sodium salt solution	1. Identification Product Name: Tristel Duo for Surfaces – Base Solution Chemical Type: Organic acid blend
2. Hazards Identification: No specific hazards Contact with acids liberates chlorine dioxide	2. Hazards Identification: No specific hazards
3. Composition/information on ingredients: Major ingredients: 0.5% Sodium Chlorite solution 2.0% Foaming Surfactant in de-mineralised water	3. Composition/information on ingredients: Major ingredients: 5% solution of Citric Acid 0.05% Boric Acid 0.05% Sorbic Acid 2.0% Foaming surfactant in de-mineralised water
4. First-aid measures: Inhalation: Non-toxic Eye contact: Rinse eyes with water Skin contact: Wash affected area with water Ingestion: Do not induce vomiting Give milk or water to drink Seek medical advice where necessary	4. First-aid measures: Inhalation: Non-toxic Eye contact: Rinse eyes with water Skin contact: Wash affected area with water Ingestion: Do not induce vomiting Give milk or water to drink Seek medical advice where necessary
5. Fire fighting measures: Non flammable	5. Fire fighting measures: Non flammable
6. Accidental release measures: Environmental Precautions: Environmental precautions required but product is biodegradable under OECD conditions operational 6/1995 Clean up method: Dispose of to clinical waste Clothing for disposal: Wear appropriate gloves	6. Accidental release measures: Environmental Precautions: Environmental precautions required but product is biodegradable under OECD conditions operational 6/1995 Clean up method: Dispose of to clinical waste Clothing for disposal: Wear appropriate gloves

This document
conforms with the
Regulations of
REACH, 2007

Tristel Duo for Surfaces Foam Foam (Activator) & Foam (Base) Safety Data Sheet

Product Name: Tristel Duo for Surfaces – Activator Solution	Product Name: Tristel Duo for Surfaces – Base Solution
7. Handling and storage:	7. Handling and storage:
<p>Handling guidelines: Safe handling in accordance with label instructions</p> <p>Not to be mixed with other chemicals</p> <p>Keep from children</p> <p>Storage guidelines: Store within a temperature range of 10° C to 35° C</p> <p>Shelf life – two years – see product for expiry date</p>	<p>Handling guidelines: Safe handling in accordance with label instructions</p> <p>Not to be mixed with other chemicals</p> <p>Keep from children</p> <p>Storage guidelines: Store within a temperature range of 10° C to 35° C</p> <p>Shelf life – two years – see product for expiry date</p>
8. Exposure controls/personal protection:	8. Exposure controls/personal protection:
<p>Personal protection: Avoid eye contact and prolonged skin contact</p> <p>Gloves and apron recommended</p> <p>Skin contact: Low risk</p> <p>Eye contact: Low risk</p> <p>Inhalation: Low risk</p> <p>Ingestion: Low risk, substantial ingestion may cause discomfort to mouth and digestive tissues</p>	<p>Personal protection: Avoid eye contact and prolonged skin contact</p> <p>Gloves and apron recommended</p> <p>Skin contact: Low risk</p> <p>Eye contact: Low risk</p> <p>Inhalation: Low risk</p> <p>Ingestion: Low risk, substantial ingestion may cause discomfort to mouth and digestive tissues</p>
9. Physical and chemical Properties:	9. Physical and chemical Properties:
<p>Physical state: Liquid</p> <p>Appearance and odour: Clear, no odour</p> <p>Evaporation rate: As water</p> <p>Boiling point: As water</p> <p>Freezing point: As water</p> <p>% Volatile (by weight): Not known</p> <p>Solubility in water (20° C): Soluble</p> <p>pH: 12.0 approximately</p> <p>Specific gravity: 1.005 @ 20° C</p>	<p>Physical state: Liquid</p> <p>Appearance and odour: Blue, no odour</p> <p>Evaporation rate: As water</p> <p>Boiling point: As water</p> <p>Freezing point: As water</p> <p>% Volatile (by weight): Not known</p> <p>Solubility in water (20° C): Soluble</p> <p>pH: 2.5 approximately</p> <p>Specific gravity: 1.025 @ 20° C</p>
	This document conforms with the Regulations of REACH, 2007

Tristel Duo for Surfaces Foam Foam (Activator) & Foam (Base) Safety Data Sheet

Product Name:	Tristel Duo for Surfaces – Activator Solution	Product Name:	Tristel Duo for Surfaces – Base Solution
15. Regulatory information:		15. Regulatory information:	
Not a licensed medicine		Not a licensed medicine	
16. Other information:		16. Other information:	
Safety phrases		Safety phrases	
(2) Keep out of reach of children		(2) Keep out of reach of children	
(24/25) Avoid contact with skin and eyes		(24/25) Avoid contact with skin and eyes	
(50) Do not mix with other chemicals		(50) Do not mix with other chemicals	
			This document conforms with the Regulations of REACH, 2007

The information and recommendations in this publication are to the best of our knowledge, information and belief accurate at date of publication. Nothing herein is to be construed as a warranty, express or implied. Final determination of the suitability of any product is the sole responsibility of the user.